

EFFICIENT AND SAFE DISINFECTION IN THREE STEPS

1 I Diffusion

Microdrop phase

Microdrop Technology®**: Biocide is sprayed in microdroplets, forming a nonwetting fog that covers all accessible and non-accessible surfaces, with excellent distribution.

2 I Contact time

Vapor phase

Evaporation of small droplets on contact with surfaces. Penetration of the biocide into the microbial cell and irreversible lethal action (full biocide spectrum).

3 I Aeration

Low H_2O_2 concentration: fast ventilation.



Why choose Phileas® Control?

Coordonated process with or without data exchange Modbus Communication RJ45 connection Industrial PLC Siemens screen

Maximum convenience to set up the disinfection process for operators Pre-start safety Disc protection sensor, notification if biocide quantity is insufficient

Biocide tanks $2 \times 2 L$ canisters, weight sensor under each canister RFID Reader biocide data (batch no., minimum durability date) Sensor H_2O_2 , temperature and humidity (Vaisala) optional

Traceability, data integrity

Audited compliant 21CFR Part 11 and GMP Annex 11 operator identification at start-up, profiles and associated rights, audit trail Diffusion reports in pdf, with compliance confirmation

Technical characteristics

Treated volume up to 300 m3*

Device flowrate 1200 mL/hr

Size (L x W x H) 650 x 800 x 1000 mm

Weight 30 kg, 4 steering wheels

Mains socket CE conformity

Minimum maintenance robust equipment

Autonomous Airborne Surface Disinfection (ASD) device, with patented technology of non-wetting fog using the centrifugation technique. In compliance with EN 17-272 standards. Marketing Authorisation AMM N° FR-2019-0071 obtained with disinfectant O2SAFE7.4®.

PERFORMANCE, SAFETY & EASE OF USE FOR THE FOLLOWING APPLICATIONS ACTIVITY SECTORS:

- Pharmaceutical industry, animal research laboratory
- Controlled access zones, aseptic production zone

^{*}according to the calculation tank capacity / 12mL/m3 dose (Marketing Authorisation)
** Microdrop Technology: micro drops resulting from the centrifugation

PERFORMANCE

- Microorganisms reduction by 4, 5 and 6 log10 depending on target (Geobacillus stearothermophilus spores)
- Effectiveness validated by Marketing Authorization N°FR-2019-0071 and in compliance with EN 17-272 standards
- Siemens PLC, touch interface and industrial environment
- Preventive maintenance contract: performance assurance for every disifnfection

SAFETY

- Operator safety: no handling of biocide, full use of canisters
- Safety of cleanroom equipment and materials: treatment with low concentration of H_2O_2 , non-corrosive, respectful of installations.
- · Security: disc protection removed before diffusion
- Security: absence or insufficient quantity of biocide before diffusion.
- Bright visibility (LEDs on dispensing modules, illuminated panels behind canisters)
- · Operator identification on the equipment
- Easy, automatic programming, reducing the risk of errors
- · Compact equipment, limiting the need for cables during disinfection cycles

EASE OF USE

- · Simple centrifugation technology, easy maintenance
- Robust, reliable equipment
- · Easy operator training
- · Equipment designed for easy cleaning, maintenance and storage
- Durability



INTEGRATION E TRACEABILITY

- Ethernet connection (RJ45 port), modbus communication with PLC and building management system data integration, communication with HVAC, doors, etc.
- H₂O₂ sensor, temperature and hygrometry (Vaisala - optional): additional cycle compliance data
- Weight sensor under each canister: integration of diffusion data

Traceability, data integrity:

- Diffusion report signed compliant, audit trail
- Equipment audited in compliance with GMP Annex 11 and 21CFR Part 11
- RFID reader: batch number and minimum durability date of biocide reported on the report, without risk of data entry error

WHY CHOOSE DEVEA?

- Biodecontamination expertise recognized in the most demanding sectors
- Range of products and services adapted to all configurations
- Market Authorization issued

Devea Biodecontamtion Services offer:

- Qualifications (IQ, OQ)
- Cycle developpement, protocol, support & performance qualification
- Training
- · Disinfection service provision
- Annual preventive maintenance contract
- Biological and chemical indicators supply
- After-sales service & equipment diagnostics within 48 hours
- Designed, manufactured and maintained in France